Double Timepix Quad + 3 GEM Module


anode plane


GEMs

readout plane

quad-boards reinforcement of anode plane

redframe












Martin Killenberg (University of Bonn)

Readout Scheme







EUDAQ: The EUDET Data Acquisition System




DataProducer: Pixelman plugin which communicated with EUDAQ

- Receives commands from Run Control
- Sends data to Data Collector
- Sends messages to Log Collector

DataCollector:

- Receives raw data
- Performs event building
- New: Plugin mechanism LCIO converter plugin for every raw data format





The Hardware Shutter



- Timepix measures clock cycles to end of shutter window
 ⇒ precise length, small jitter
- Software shutter is not precise (Windows timing: 1ms)
 ⇒ use precise hardware pulse generator

Software problem:


- Latest Pixelman version does not support hardware shutter any more
 Bypass the Muros and give shutter signal directly to chip
- Attention: Shutter and readout have to be synchronised. Opening the shutter while reading out the chip will disrupt the communication. Chip has to be reset.
 - TLU inhibits triggers while device is busy.





The TLU Trigger-Busy-Handshake





- TLU signalises trigger for each device
- DAQ device asserts busy, TLU takes down trigger
- TLU blocks all further triggers while any of the devices is busy



Busy 2

TLU Problem



The TLU does not always block the trigger properly:



- Event numbers get out of sync
- Communication with Timepix chip is disrupted




Solution to the TLU Problem



Special "OR" realised using a programmable logic box (NEMBox):

- Channels are raised asynchronously
- Channels are lowered synchronously



 $\Rightarrow \mathsf{Trigger} \; \mathsf{system} \; \mathsf{works} \; \mathsf{perfectly!}$





Ar/CO_2 70/30 (finished):

- Different drift distances
- Different angles
- Different particle energies
- Different GEM settings

T2K Gas (starts tomorrow):

- Different drift distances
- Different angles
- Different particle energies
- Different GEM settings
- Laser dots



