

Direct Coupling of SiPMs to Scintillator Tiles: Timing and Uniformity

Christian Soldner

Max-Planck-Institut für Physik

Tungsten Prototype Workshop
DESY, 2. März 2010



Max-Planck-Institut für Physik
(Werner-Heisenberg-Institut)





Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary

- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand
- 3 Timing of the Tile Signal
- 4 Determination of the Tile (Non-)Uniformity
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary

- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand
- 3 Timing of the Tile Signal
- 4 Determination of the Tile (Non-)Uniformity
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



The CALICE Analog Hadron Calorimeter (AHCAL)

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity



Tile
Modification

Summary

Development of blue sensitive SiPMs

allow for photon readout without wavelength shifting fiber:

- **Advantage:** Avoid machining of WLS into millions of tiles (ILD)
Avoid difficult alignment of SiPM to WLS-end
Improve timing of light collection
- **Non-uniform cell readout:** Measured energy deposition varies with the lateral position at which a particle traverses the tile
- Position dependence affects the performance of the HCAL.





The Effect of Non-Uniform Cell Readout

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary

Consequences: Subject of further investigations

Consequence 1

Distortion of the cell signal distribution
⇒ Effect on HCAL calibration

Consequence 2

Distortion of the energy sum distribution of showers
⇒ Effect on energy reconstruction process

Consequence 3

Deterioration of the energy resolution increasing with the particle energy



Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

**The Test
Stand**

Timing

Uniformity

Tile
Modification

Summary

- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand**
- 3 Timing of the Tile Signal
- 4 Determination of the Tile (Non-)Uniformity
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



The Experimental Setup

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation


The Test
Stand

Timing

Uniformity

Tile
Modification

Summary



Test Stand Properties

- Radioactive ^{90}Sr source: Beta decay ($E_{\text{end point}} = 2.27 \text{ MeV}$)
- Movable stage: Translation in XY-direction over $3 \times 3 \text{ cm}^2$ tile
- Active (Fan, Air Cond.) and Passive (Thermistor) T-control
- Select tile traversing electrons \rightarrow Active coincidence trigger



Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary



- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand
- 3 Timing of the Tile Signal**
- 4 Determination of the Tile (Non-)Uniformity
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



Time Resolved Measurements

Fast timing of HCAL implies fast light collection within tiles

- Acquire waveforms of SiPM Signal with 4 GHz Oscilloscope
- High sampling: Arrival of every single photon on SiPM can be identified



Direct Coupling

Signal from directly coupled tile significantly faster:
No delay due to absorption and reemission in WLS fiber



Time Resolved Measurements

SiPMs and Scintillator Tiles

C. Soldner

Motivation



The Test Stand

Timing

Uniformity

Tile Modification

Summary



CALICE 1st generation tile with WLS fiber

- Broad signal peak
- Long integration times needed

Directly Coupled Tile

- Fast peaking signal, pronounced peak
- Short integration times sufficient



Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary

- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand
- 3 Timing of the Tile Signal
- 4 Determination of the Tile (Non-)Uniformity**
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



The Lateral Uniformity Scan

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing


Uniformity

Tile
Modification

Summary

Scanning Sequence

- First test: Direct SiPM coupling to one side of a $3 \times 3 \text{ cm}^2$ tile
- Recording of $500 e^-$ signals at 60×60 XY-positions on the tile
→ Pixel Size: $0.5 \times 0.5 \text{ mm}^2$
- Determine mean signal height (MSH) from signal distribution
- MSH vs. XY-Position: Observe strong non-uniformity
- Modification of the tile geometry necessary to restore uniformity





The Lateral Uniformity Scan

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing



Uniformity

Tile
Modification

Summary

Scanning Sequence

- First test: Direct SiPM coupling to one side of a $3 \times 3 \text{ cm}^2$ tile
- Recording of 500 e^- signals at 60×60 XY-positions on the tile
→ Pixel Size: $0.5 \times 0.5 \text{ mm}^2$
- Determine mean signal height (MSH) from signal distribution
- MSH vs. XY-Position: Observe strong non-uniformity
- Modification of the tile geometry necessary to restore uniformity





The Lateral Uniformity Scan

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing



Uniformity

Tile
Modification

Summary

Scanning Sequence

- First test: Direct SiPM coupling to one side of a $3 \times 3 \text{ cm}^2$ tile
- Recording of $500 e^-$ signals at 60×60 XY-positions on the tile
→ Pixel Size: $0.5 \times 0.5 \text{ mm}^2$
- Determine mean signal height (MSH) from signal distribution
- MSH vs. XY-Position: Observe strong non-uniformity
- Modification of the tile geometry necessary to restore uniformity





The Lateral Uniformity Scan

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing



Uniformity

Tile
Modification

Summary

Scanning Sequence

- First test: Direct SiPM coupling to one side of a $3 \times 3 \text{ cm}^2$ tile
- Recording of 500 e^- signals at 60×60 XY-positions on the tile
→ Pixel Size: $0.5 \times 0.5 \text{ mm}^2$
- Determine mean signal height (MSH) from signal distribution
- MSH vs. XY-Position: Observe strong non-uniformity
- Modification of the tile geometry necessary to restore uniformity





The Lateral Uniformity Scan

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing



Uniformity

Tile
Modification

Summary

Scanning Sequence

- First test: Direct SiPM coupling to one side of a $3 \times 3 \text{ cm}^2$ tile
- Recording of 500 e^- signals at 60×60 XY-positions on the tile
→ Pixel Size: $0.5 \times 0.5 \text{ mm}^2$
- Determine mean signal height (MSH) from signal distribution
- MSH vs. XY-Position: Observe strong non-uniformity
- **Modification of the tile geometry necessary to restore uniformity**





Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

**Tile
Modification**

Summary

- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand
- 3 Timing of the Tile Signal
- 4 Determination of the Tile (Non-)Uniformity
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



Towards an Optimal SiPM-Tile Entity

SiPMs and Scintillator Tiles

C. Soldner

Motivation



The Test Stand

Timing

Uniformity

Tile Modification

Summary



Tile Development

Simple Tile: SiPM coupled to the side

Highly non-uniform

⇒ Solution: Reduce scintillating material close to SiPM

Overall mean of signal height: 13.0 p.e.

⇒ Solution: Integrate SiPM into the tile



Towards an Optimal SiPM-Tile Entity

SiPMs and Scintillator Tiles

C. Soldner

Motivation



The Test Stand

Timing

Uniformity

Tile Modification

Summary



Tile Development

Integrated Dimple: Height 1 mm, Width 4 mm, Depth 5.5 mm

⇒ Overall mean: 18.4 p.e. (Diffuse light refraction at dimple)

⇒ Good uniformity

⇒ Signal at SiPM position not lower than 13 p.e.

⇒ Scalable solution for mass production

(arXiv:1001.4665 [physics.ins-det])



Towards an Optimal SiPM-Tile Entity

SiPMs and Scintillator Tiles

C. Soldner

Motivation



The Test Stand

Timing

Uniformity

Tile Modification

Summary



Tile Development

Transferable concept for the next generation prototype: 3mm Option:

⇒ Overall mean: 13.15 p.e.

⇒ Signal at Dimple not lower than 8 p.e.



Towards an Optimal SiPM-Tile Entity

SiPMs and Scintillator Tiles

C. Soldner

Motivation



The Test Stand

Timing

Uniformity

Tile Modification

Summary



Tile Development

Very small packaging of new MPPCs

→ opens up the possibility for new tile geometries:

Aim: Avoid Signal Drop at SiPM Coupling Position

→ Very small integration hole

→ 2 mm deep bottom dimple (spherical drilling head with 5 mm radius)

Note: Only tested for 5 mm tiles, SiPM is SMD



High Precision Uniformity: Side Coupling vs. Integrated Dimple Concept

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity


Tile
Modification

Summary

Side Coupling

Part of the Scintillator Tile	Deviation of overall mean: 13.0
98% (91%)	$\pm 20\%$
94% (81%)	$\pm 10\%$
69% (57%)	$\pm 5\%$

Unprecise → Cut Tile Edges!





High Precision Uniformity: Side Coupling vs. Integrated Dimple Concept

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary



Side Coupling

Part of the Scintillator Tile	Deviation of overall mean:
98% (91%)	$\pm 20\%$
94% (81%)	$\pm 10\%$
69% (57%)	$\pm 5\%$

Unprecise \rightarrow Cut Tile Edges!

Side Dimple

Part of the Scintillator Tile	Deviation of overall mean:
99% (90%)	$\pm 20\%$
97% (84%)	$\pm 10\%$
88% (73%)	$\pm 5\%$





High Precision Uniformity: Side Coupling vs. Integrated Dimple Concept

SiPMs and
Scintillator
Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary



Side Coupling

Part of the Scintillator Tile	Deviation of overall mean:
98% (91%)	$\pm 20\%$
94% (81%)	$\pm 10\%$
69% (57%)	$\pm 5\%$

Unprecise → Cut Tile Edges!

Bottom Dimple

Part of the Scintillator Tile	Deviation of overall mean:
99.8% (90%)	$\pm 15\%$
96% (82%)	$\pm 10\%$
83% (69%)	$\pm 5\%$





Outline

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test
Stand

Timing

Uniformity

Tile
Modification

Summary

- 1 Motivation: Reasons for Direct Coupling
- 2 The Experimental Test Stand
- 3 Timing of the Tile Signal
- 4 Determination of the Tile (Non-)Uniformity
- 5 Optimization of the Scintillator Tile Geometry
- 6 Summary and Outlook



Summary and Outlook

SiPMs and Scintillator Tiles

C. Soldner

Motivation

The Test Stand

Timing

Uniformity

Tile Modification

Summary

Timing of the Tile Signal

Observe directly: Significantly faster tile signal for direct coupling

Tile Uniformity Study

- SiPM Integration: Tight cell mounting, higher light yield, easy SiPM alignment
- Side Dimple: High Uniformity, easy design, scalable for 5 mm and 3 mm tiles
- NEW: Bottom Dimple: Only possible with very small casing
→ MPPC P-Series available as SMD type
Avoids signal drop at coupling pos while keeping light yield high
→ Potential for perfect uniformity
Further improvements under investigation!
Test concept for 3 mm tiles!



Outline

SiPMs and
Scintillator
Tiles

C. Soldner

Appendix

7 Appendix





The Experimental Setup

SiPMs and
Scintillator
Tiles

C. Soldner

Appendix



Trigger condition

Signal of Tile & Coincidence Cube simultaneously > 3 p.e.



Tile Coupling




SiPMs and
Scintillator
Tiles

C. Soldner

Appendix

Coupling Properties

SiPM:	MPPC25, 1600 Pix, 2 plastic casings
Scintillator:	Saint Gobain BC-420
Tile Size:	$3 \times 3 \times 0.5$ and $3 \times 3 \times 0.3$ cm ³
Tile Surface:	Polished, completely enclosed by 3M mirror foil
SiPM Coupling:	Direct (Air Gap) fiberless coupling
Coupling Pos:	Center of one side face






Quantification of the Uniformity with Area Fractions

SiPMs and
Scintillator
Tiles

C. Soldner

Appendix



Determination of the Area Fractions

- Define the Tile Position (fixed to $3 \times 3 \text{ cm}^2$):
→ Note: Measurement unprecise at tile edges
Reason: Finite iris size of source casing → conelike e^- emission
- Determine OMSH: Mean value of all measurement points outside the extreme regions
- Determine Area Fractions: Fraction of all on-tile measurement points lying within a certain deviation region around the OMSH




Quantification of the Uniformity with Area Fractions

SiPMs and
Scintillator
Tiles

C. Soldner

Appendix



Determination of the Area Fractions

- Define the Tile Position (fixed to $3 \times 3 \text{ cm}^2$):
→ Note: Measurement unprecise at tile edges
Reason: Finite iris size of source casing → conelike e^- emission
- Determine OMSH: Mean value of all measurement points outside the extreme regions
- Determine Area Fractions: Fraction of all on-tile measurement points lying within a certain deviation region around the OMSH




Quantification of the Uniformity with Area Fractions

SiPMs and
Scintillator
Tiles

C. Soldner

Appendix



Determination of the Area Fractions

- **Define the Tile Position (fixed to $3 \times 3 \text{ cm}^2$):**
→ Note: Measurement unprecise at tile edges
Reason: Finite iris size of source casing → conelike e^- emission
- **Determine OMSH:** Mean value of all measurement points outside the extreme regions
- **Determine Area Fractions:** Fraction of all on-tile measurement points lying within a certain deviation region around the OMSH



The Effect of Non-Uniform Cell Readout

SiPMs and
Scintillator
Tiles

C. Soldner

Appendix

Consequences: Needs further investigation

Consequence 1

Distortion of the cell signal distribution
⇒ Effect on HCAL calibration

Consequence 2

Distortion of the energy sum distribution of showers
⇒ Effect on energy reconstruction process

Consequence 3

Deterioration of the energy resolution increasing with the particle energy